





Catalogue # 3479 Pancreas Carcinoma of Head

Description	Image
Fluid dark coffee color. Cannot see specimen. Fluid unknown.	
When specimen was removed from container residual mold was detected encrusting the specimen. ARROWS Fluid had a phenol smell.	
	

Catalogue # 3479 Pancreas Carcinoma of Head

Description	Image
<p>Support glass rod for mounting broken. Attempted fabrication of new mounting frame from acrylic. Acrylic was stable but glue joints failed. A glass mounting plate was cut and drilled for attachment sutures</p>	
<p>Rinsed specimen with DI water, soaking for 2x changes. 10/25 → fresh DI water 10/26 → 15% ETOH 11/1 → 30% ETOH 11/18 → 50% ETOH 11/1 → 70% ETOH</p> <p>Cleaned container, restored old label, remounted specimen, resealed lid.</p>	